


TheraCem[®]

CEMENTATION. THE THERA WAY.



Rx Only

TheraCem[®]

Self-Adhesive Resin Cement

The NEXT GENERATION of Cements!

TheraCem is a dual-cured, calcium and fluoride-releasing, self-adhesive resin cement indicated for luting crowns, bridges, inlays, onlays and posts (prefabricated metal/nonmetal/fiber posts). It is the first BISCO cement to incorporate the THERA technology.

TheraCem has continuous release of calcium and fluoride ions.¹ It is at first acidic, which is needed for initial bonding, but becomes alkaline within minutes.²

The MDP monomer containing formula promotes adhesion that enhances the bond strength to zirconia, metal, and the tooth structure without the use of an additional dental adhesive or primer. Additionally, it offers a high degree of conversion for increased strength and high radiopacity for differentiation on radiographs.

TheraCem's single step procedure makes cementation to zirconia easier than ever before.



Calcium & Fluoride Releasing
Continuous ion release¹



High Degree of Conversion
Ensures a higher physical strength



Alkaline pH
Transitions from acidic to an alkaline pH in minutes²



Easy to Use
Auto-mix, dual-syringe provides a consistent mix for immediate delivery



Easy Clean-Up
Specifically formulated to allow for quick and easy removal of cement



Strong Bond to Zirconia
Delivers a strong bond to zirconia and most substrates with no priming or etching required



MDP
TheraCem contains MDP which allows to bond to zirconia, metal and alumina substrates without the use of an additional primer



Dual-Cured
Material will cure even in deep restorations where light cannot reach



Self-Adhesive
No bonding agents required



Highly Radiopaque
Easy to identify on radiographs for quick and effective diagnosis



Data on file. BISCO, Inc.

1. Gleave CM, Chen L, Suh BI. Calcium & fluoride recharge of resin cements. Dent Mater. 2016 (32S):e26.

2. New Self-adhesive Resin Cement With Alkaline pH. Chen L, Gleave C, Suh B, J Dent Res 96(A):#286, 2017.

■ DUAL-CURED

■ EASY APPLICATION

■ RADIOPAQUE

■ ALKALINE pH

THERA HYDROPHILIC MATRIX




BISCO has developed a hydrophilic matrix that allows for ion exchange. Previous or traditional resin matrices have been hydrophobic, but BISCO's matrix allows for ion exchange (as water goes into the matrix, whereby calcium hydroxide ions and fluoride ions are released).


TheraCem Technique

Glass Ceramics / Porcelain (silica-based)
or Lithium Disilicate / Silicate


Zirconia, Metal,
Alumina, Composite



TheraCem Continuous Fluoride and Calcium Release¹



Shear Bond Strength to Zirconia



1. Gleave CM, Chen L, Suh BI. Calcium & fluoride recharge of resin cements. Dent Mater. 2016 (32S):e26.

TheraCem Case

Dentistry courtesy of Dr. Jordan Sebastian DDS and Dr. Burgos Hugo DDS




Cast metal post and core cemented on the first upper left bicuspid.




A CAD/CAM zirconia restoration was fabricated.



Restoration tried-in to check proper fit and occlusion.



ZirClean® applied for 20 seconds to remove surface contaminants, then rinsed and dried.



TheraCem dispensed into the intaglio surface of the restoration.



The restoration was seated, tack cured, and excess cement removed.



Final restoration.

Ordering Information

TheraCem Single Dual-Syringe Package

1 Dual-Syringe (8g), 15 Mixing Tips and Intraoral Tips, Instructions

Natural Shade.....D-46311P

Exclusively distributed by Curion

curion.ca | EN: 1.800.667.8811 | FR: 1.800.211.1200

 **CURION**
formerly Bisco Canada