



# Smart Fitness Scale

**Model:** FIT 8S

## User Manual



### Questions or Concerns?

Mon–Fri, 9:00 am–5:00 pm PST/PDT  
support@etekcity.com • (855) 686-3835

# Thank you for purchasing the Smart Fitness Scale by Etekcitey.

If you have any questions or concerns, please reach out to our helpful Customer Support Team at **[support@etekcity.com](mailto:support@etekcity.com)**.  
We hope you enjoy your new smart fitness scale!

# Table of Contents

|                       |           |
|-----------------------|-----------|
| Package Contents      | <b>4</b>  |
| Specifications        | <b>4</b>  |
| Safety Information    | <b>5</b>  |
| Function Diagram      | <b>8</b>  |
| Digital Display       | <b>9</b>  |
| VeSync App Setup      | <b>10</b> |
| Before First Use      | <b>10</b> |
| Using the Smart Scale | <b>11</b> |
| Using the App         | <b>12</b> |
| Maintenance           | <b>13</b> |
| Troubleshooting       | <b>16</b> |
| Warranty Information  | <b>18</b> |
| Customer Support      | <b>21</b> |

# Package Contents

- 1 x Smart Fitness Scale
- 3 x 1.5V AAA Batteries (Pre-Installed)
- 1 x Quick Start Guide

# Specifications

|                        |  |
|------------------------|--|
| <b>Weight Capacity</b> | 400 lb / 28 st / 180 kg                    |
| <b>Increments</b>      | 0.1 lb / 0.05 kg                           |
| <b>Weight Units</b>    | lb / st / kg                               |
| <b>Platform</b>        | Tempered glass                             |
| <b>Dimensions</b>      | 10.2 x 10.2 x 0.8 in /<br>26 x 26 x 2.1 cm |
| <b>Battery</b>         | 3 x 1.5V AAA Batteries                     |
| <b>Auto-Off</b>        | 30 seconds                                 |



# READ AND SAVE THESE INSTRUCTIONS

## Safety Information

**Please read and follow all instructions and safety guidelines in this manual.**

### Medical

- This scale should not be used to diagnose or treat any medical condition.
- Scale data should not be used as a substitute for medical advice from a healthcare professional.
- This scale uses a harmless electrical current to measure body fat. Consult a healthcare professional before using the smart scale if you use medical electronic equipment such as a pacemaker.

### General Use

- Check the scale before using. **Do not** use the scale if it is damaged in any way.
- **Do not** strike or hit the scale. **Do not** use the scale if the glass is cracked or broken. Due to the nature of tempered glass, the glass will shatter if there is even a slight crack. If you see a crack, contact **Customer Support** immediately (page 21).



- The glass platform is slippery when wet. **Always** make sure the platform and your feet are dry before stepping on for measurement.
- Place the scale on a hard, flat surface to avoid tipping.
- To weigh safely and accurately, stand with your feet shoulder-width apart while on the scale. **Do not** stand on the edge of the scale or jump on the scale.
- Avoid standing on the scale for a long period of time, or when not weighing yourself.
- **Do not** overload the scale. Maximum weight capacity is 400 lb / 28 st / 180 kg.
- Keep the scale in a cool area with low humidity. Keep the scale away from heat sources, such as ovens or heaters. Avoid temperature fluctuations. Store the scale in a horizontal position, not a vertical position.
- Handle with care.
- Not for commercial use. Household use **only**.




## Batteries

- **Do not** mix new and old batteries.
- If the low battery indicator (“Lo”) appears on the display, replace the batteries.
- If the scale is not in use for a long period of time, remove the batteries to extend the lifetime of the scale.



# Function Diagram



**A.** Weighing Platform

**B.** Display

**C.** Anti-Skid Padding


**D.** Unit Switching

Button

**E.** Battery Compartment



# Digital Display



- A.** Weight Measurement
- B.** Wireless Connection Icon
- C.** Weight Units



# VeSync App Setup

**Note:** Due to app updates and improvements, the contents of this manual may change without notice. The VeSync app user interface may appear slightly different.

1. To download the VeSync app, scan the QR code or search “VeSync” in the Apple App Store® or Google Play™ Store.

**Note:**

- This scale uses the VeSync app, not the VeSyncFit app.
  - For Android™ users, you may need to select **Allow** to use VeSync.
2. Open the VeSync app. If you already have an account, tap **Log In**. To create a new account, tap **Sign Up**.



## Before First Use

1. Open the battery compartment on the underside of the scale. Remove the included AAA batteries.
2. Remove the plastic insulator strip from the battery compartment.
3. Replace the batteries. Follow the guidelines in the battery compartment to match the **+** and **-** terminals.




# Using the Smart Scale

## Note:


- To switch between kilograms (kg), stone (st), and pounds (lb), press the unit switching button on the underside of the scale.
- If you are pregnant, scale measurements will not be accurate.
- If you want to measure the weight of your baby or pet, you can click the “<” (Expand) button in the lower right corner of the home page and click the “😊” (Baby Mode) button to enter the baby measurement mode and measure according to the prompts. In baby mode, the APP will only measure weight and BMI.

## Weighing Without the VeSync App

1. Place the scale on a hard, flat surface. [\[Figure 1.1\]](#)
2. Step on the scale, then step off. Wait for the screen to display “0.0” before using.
3. Place dry, bare feet on the scale’s conductive areas. Stand on the scale until the display flashes several times, showing your final measurement. [\[Figure 1.2\]](#)



[\[Figure 1.1\]](#)



[\[Figure 1.2\]](#)



# Using the App

## Creating a User

If you are a new user, you will need to enter your physical information before weighing yourself for the first time so the scale can calculate your body metrics.

**Note:**

- *Entering inaccurate information will make the scale's measurements less accurate.*
- *If you don't identify with the available gender options, or are unsure, pick the choice that you think best physically matches you. This choice is used for physical measurements, such as body fat percentages.*




# Maintenance

## Clearing the Scale

If the scale has been moved or flipped upside down, it must be cleared to ensure accurate results.

1. Place the scale on a hard, flat surface.
2. Step on the scale until digits appear on the display, then step off. [\[Figure 2.1\]](#)
3. The scale will show “0.00” when cleared. [\[Figure 2.2\]](#)



[\[Figure 2.1\]](#)



[\[Figure 2.2\]](#)



## Cleaning the Scale

- Use a slightly damp cloth to clean the surface of the scale, followed by a soft, dry cloth to wipe away any remaining moisture.
- **Do not** use chemical or abrasive cleaning agents to clean the scale.

## Replacing the Batteries

The display will show “**Lo**” when the batteries need to be replaced.

1. Open the battery compartment on the underside of the scale.
2. Remove the old batteries and dispose of them properly.
3. Install 3 new 1.5V AAA batteries. Follow the guidelines in the battery compartment to match the **+** and **-** terminals.
4. After closing the battery compartment, clear the scale (see **Clearing the Scale**, page 13).



# Display Messages

A 7-segment display showing the number '68.0'. The first digit is '6', the second is '8', and the third is '0'. A small 'e' symbol is positioned above the top-right corner of the '0'.

Wireless Connection

A 7-segment display showing the characters 'Lo'. The 'L' is formed by the left vertical bar, the bottom horizontal bar, and the bottom-right vertical bar. The 'o' is formed by the top horizontal bar, the right vertical bar, and the bottom horizontal bar.

Change Battery

A 7-segment display showing the character 'E'. It is formed by the left vertical bar, the top horizontal bar, and the bottom horizontal bar.

Clearing

A 7-segment display showing the characters 'Err'. The 'E' is formed by the left vertical bar, the top horizontal bar, and the bottom horizontal bar. The first 'r' is formed by the right vertical bar, the top horizontal bar, and the bottom horizontal bar. The second 'r' is formed by the right vertical bar and the bottom horizontal bar.

Overload



# Troubleshooting

## **The scale doesn't turn on.**

- Check if the batteries are properly installed. If using the scale for the first time, remove the plastic insulator strip from the battery compartment.

## **Information flashes on the screen and then disappears.**

- The batteries are low and may need to be replaced.

## **The scale doesn't connect with the VeSync app.**

- Make sure your phone's operating system is running iOS® 10.0+ or Android™ 5.0+.
- Close and relaunch the VeSync app. Make sure you're using the latest version.
- Turn on Bluetooth® on your phone (for Android users, make sure Location is turned on as well).
- Make sure your scale and phone are within 10 m / 30 ft of each other.
- Make sure the scale is not currently connected to any other phone or app.



## **My scale only measures my weight, not my body fat or other metrics.**

- Set up your smart scale in the VeSync app. Your weight will show on the scale display, and more results will show on the app.
- Make sure Bluetooth is turned on.
- If you're wearing shoes or socks, the scale cannot measure health metrics such as body fat. Step on the scale with bare feet.

**Note:** *If the scale has errors on the display, will not turn off automatically, or has a problem not listed here, try removing the batteries for 3 seconds, then replacing them.*

If you still need help, please contact **Customer Support** (page 21).



# Federal Communication Commission Interference Statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) This device must accept any interference received, including interference that may cause undesired operation.

**NOTE:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

**FCC Caution:** Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.

## FCC Radiation Exposure Statement

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. End users must follow the specific operating instructions for satisfying RF exposure compliance. To maintain compliance with FCC RF exposure compliance requirements, please follow operation instructions as documented in this manual. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter. The availability of some specific channels and/or operational frequency bands are country dependent and are firmware programmed at the factory to match the intended destination. The firmware setting is not accessible by the end user.



# Warranty Information

|   |                     |
|---|---------------------|
| <b>Product</b>  | Smart Fitness Scale |
| <b>Model</b>  | FIT 8S              |
| For your own reference, we strongly recommend that you record your order ID and date of purchase. |                     |

## TERMS & POLICY

Etekcitey Corporation (“Etekcitey”) warrants this product to the original purchaser to be free from defects in material and workmanship, under normal use and conditions, for a period of **2 years** from the date of original purchase.

Etekcitey agrees, at our option during the warranty period, to repair any defect in material or workmanship or furnish an equal product in exchange without charge, subject to verification of the defect or malfunction and proof of the date of purchase.

There is no other express warranty. This warranty does not apply:

- If the product has been modified from its original condition;
- If the product has not been used in accordance with directions and instructions in the user manual;
- To damages or defects caused by accident, abuse, misuse or improper or inadequate maintenance;
- To damages or defects caused by service or repair of the product performed by an unauthorized service provider or by anyone other than Etekcitey;
- To damages or defects occurring during commercial use, rental use, or any use for which the product is not intended;
- To damages or defects exceeding the cost of the product.

Etekcitey will not be liable for indirect, incidental, or consequential damages in connection with the use of the product covered by this warranty.

This warranty extends only to the original consumer purchaser of the product and is not transferable to any subsequent owner of the product regardless of whether the product is transferred during the specified term of the warranty.

This warranty does not extend to products purchased from unauthorized sellers. Etekcitey's warranty extends only to products purchased from authorized sellers that are subject to Etekcitey's quality controls and have agreed to follow its quality controls.

All implied warranties are limited to the period of this limited warranty.

This warranty gives you specific legal rights and you may also have other rights which vary from state to state. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

If you discover that your product is defective within the specified warranty period, please contact Customer Support via **support@etekcity.com**. **DO NOT** dispose of your product before contacting us. Once our Customer Support Team has approved your request, please return the product with a copy of the invoice and order ID.

Every Etekcitey product automatically includes a 2-year warranty. To make the customer support process quick and easy, register your product online at [www.etekcity.com/warranty](http://www.etekcity.com/warranty).

This warranty is made by:

Etekcitey Corporation  
1202 N. Miller St., Suite A  
Anaheim, CA 92806



# Customer Support

If you have any questions or concerns about your new product, please contact our helpful Customer Support Team.

## **Etekcitey Corporation**

1202 N. Miller St., Suite A  
Anaheim, CA 92806

**Email:** [support@etekcity.com](mailto:support@etekcity.com)

**Toll-Free:** (855) 686-3835

## **Support Hours**

Mon–Fri, 9:00 am–5:00 pm PST/PDT

*\*Please have your order invoice and order ID ready before contacting Customer Support.*



# Attributions

Google, Android, and Google Play are trademarks of Google LLC.

App Store® is a trademark of Apple Inc.

iOS is a registered trademark of Cisco Systems, Inc. and/or its affiliates in the United States and certain other countries.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Etekcitey Corporation is under license. Other trademarks and trade names are those of their respective owners.





**Connect with us @Etekcinity**

