

medipal

# Detergent Wipes

Ideal for multi-surface cleaning and damp dusting, Medipal Detergent Wipes have an innovative low-smear formula to provide effective cleaning with no visible residue. Alcohol free and easy-to-use.



## At a glance:




- Alcohol free
- Neutral pH
- Low smear formula
- No visible residue

## Ideal for:



- Cleaning and damp dusting non porous surfaces
- Cleaning patient shared equipment where there is a low risk of infection
- Removing soiling from surfaces before disinfection

## Material compatibility:



- Acrylic (PMMA)
- Polycarbonate (PC)
- Polythene (PE)
- Polyvinylchloride (PVC)
- Waterproof fabric
- 316 Stainless Steel
- 304S15 Stainless Steel

## Sustainability:



- Lighter 2 litre canister with 27% plastic reduction
- Tritex material, canister, buckets, lids, and labels are fully recyclable
- Cardboard used in our outer packaging is FSC approved




# medipal

## How to use


### STEP 1 Risk assessment

Please follow your agreed risk assessment policy guidelines regarding the use of PPE.




### STEP 2 Wipe selection

Choose the dispenser type and dispense the wipe.




### STEP 3 S-Shape technique

Wipe the surface in an S-Shape moving from clean to dirty. Use the wipe flat not scrunched. Do not go over the same area twice with the same wipe. Use a fresh wipe if your wipe becomes soiled or dry.




### STEP 4 Discard

Discard used wipes in the appropriate waste bin following your local agreed guidelines.



### STEP 5 Let dry

Allow the surface to dry naturally before use.



# medipal

## Technical details

Medipal's proprietary Tritex® material is made by bonding layers of spunbond and meltblown polypropylene fibres together in a high heat process. The addition of a hydrophilic coating means optimum and consistent liquid discharge.

### Chemical formulation:

Medipal Detergent cleaning solution is formulated using environmentally friendly and VOC free ingredients. The specially designed formula contains surfactants and a naturally derived ingredient which removes contamination from surfaces leaving them clean and streak free. The formulation is maintained at a skin balanced pH and is dermatologically tested to ensure skin safety for frequent users.

### Materials:

Tritex® material is made by bonding layers of spunbond and meltblown polypropylene fibres together in a high heat process before a hydrophilic coating is applied. This material is 100% synthetic.



#### High strength

The material offers high strength in both directions – across and along the wipe.



#### Excellent solution retention

Ensuring even wetness throughout the wipe.



#### Very low linting

Low risk of leaving contaminating fibres on surfaces.



#### Efficient solution release

Active wipe ingredients are transferred to the surface rather than being trapped in the wipe material fibres.

# medipal

## Products within this range



### 275 Wipe Bucket

Packs per case – 4

Product Code – W532110MP



### 350 Wipe Bucket

Packs per case – 4

Product Code – W535110MP



### 150 Wipe Canister

Packs per case – 10

Product Code – W530110MP