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FR

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EN

OPERATING INSTRUCTIONS FOR DIGITAL THERMOMETER

IMPORTANT. READ INSTRUCTIONS BEFORE USE AND KEEP FOR FUTURE REFERENCE

FEATURES	
- Temperature Range: 32 °C ~ 42,9 °C (90 °F ~ 109,9 °F).	- Battery life: about 250 hours or 1000 times.
- Accuracy : ± 0.1 °C (35 °C ~ 42 °C) ; ± 0.2 °C (< 35 °C).	- Flexible tip for comfort and safety.
- Digital LCD display.	- Beeps when reaches highest temperature.
- Battery: 1.5 VDC (LR41) alkaline.	- Auto-off in approximately 10 minutes.
- Power consumption: 0.15 mW.	- Last measurement memory.
- Dimensions: 129 mm (L) x 21 mm (W) x 12.5 mm (H).	- Low battery indicator "■".
	- Weight with battery: 12 g.

HOW TO TAKE TEMPERATURE?

- Press power switch, the screen shows the full segment 188.8 °F with a sound of "B", which indicates sufficient power and normal self-test routine as well as for measurement.
- Memory function: after the full segment, the last measured value with the letter "M" are displayed, then to show 37 °C / 98,6 °F.
- Release the button, then displays « Lo °C » or « Lo °F » and ready to measure. During measurement the symbol « °C » or « °F » will flash while the temperature is climbing. This will also happen when temperature increase remains below 0,1 °C / 0,2 °F for about 16 seconds. If temperature remains below 32 °C / 90 °F, the symbol « °C » or « °F » will flash for about 60 seconds. Then it will auto-off for saving power.
- A buzzing tone will be heard for about 10 seconds which means the finish of measurement. Normal alarm: sound of « bi...bi...bi...bi » means normal temperature that below 37,8 °C / 100 °F. Fever alarm: sound of « bi,bi,bi » means fever that the temperature is above 37,8 °C / 100 °F.

DIFFERENT METHODS OF USE

- You are advised to measure your temperature at specific time in morning, if it is earlier or later, the temperature would be lower or higher.
- **Used in mouth.** Please close your mouth for 1 minute for stability of temperature, breath calmly through nose. Respiratory air would affect the accuracy of measure of your mouth. Normal temperature: about 37 °C / 98,6 °F.
- **Used in anus of rectum.** This is the way that can get the most steady temperature and suits for both infants and children. The sensor shall be inserted for 2-3 cm. Vaseline can be used for easier inserting. Normal temperature: about 37,6 °C / 99,7 °F.
- **Used under armpit.** Clean the thermometer with a dry cloth and put it under armpit. Press and fix it with the arm. The temperature value of this measurement is not so accurate that no recommended for accurate measurement. Normal temperature: about 36,7 °C / 98,1 °F.
- **Notes:** the normal temperatures are different from the different human bodies and measuring time. Sterilization is required before every measurement.

BATTERY REPLACEMENT

Battery must be replaced if "■" symbol appears during measurement. Make sure a battery of the same type is on hand and proceed as follows: Remove battery compartment lid and change with a button battery (LR41). Ensure battery is correctly positioned with "+" facing up.

WARNING: keep battery out of children's reach, swallowing battery could be fatal, battery should not be charged or placed into extreme heat as it may explode.

CLEANING AND MAINTENANCE

The thermometer is waterproof and can be cleaned by water or sterilization lotion with soft cloth: Sterilization lotion below is recommended:
Name: Density - **Medical alcohol:** 95 %
Notes: other lotion, density or cleaning method may damage the thermometer.

CAUTIONS

- Store the unit in the protective case while not in use.
- Do not store or keep in direct sunlight, or in high temperatures.
- Keep away from strong electromagnetic source to avoid interference.
- Do not open the thermometer except to replace battery.
- Do not drop the unit, as this may ruin the unit.
- It shall be operated in an environment at the temperature of 10 °C / 50 °F and the humidity of 40-80 % while stored and transported in an environment at the temperature of -10 °C ~ 40 °C (14 °F ~ 104 °F) and humidity of 30 ~ 90 %.
- Never put it in boiled water for sterilizing.
- Please take out the battery if not use the unit for longtime.

REMARKS

- It is suitable for measuring temperature of human body and necessary of families.
- It is designed and manufactured complying with relevant requirements in EN12470-3.
- Unwanted electrical and electronic items are collected separately. Do not discard unwanted electrical and electronic items with unsorted municipal waste, but ensure that they are collected separately.

Complies with safety requirements.

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