

ZirClean™

Restoration Cleaner

ZirClean is a cleaning gel designed for the non-abrasive cleaning of the bonding surfaces of zirconia (and other prosthetic restorations) after intraoral try-in. ZirClean helps to achieve reliable adhesive cementation results by removing the phosphate contamination of zirconia (as well as ceramic and metal restoration surfaces) that occur during try-in.*

Product Highlights


- ▶ Delivery system offers improved ease-of-use
- ▶ Easy placement and clean up
- ▶ Helps achieve reliable adhesive cementation results*
- ▶ Contaminants such as saliva can be removed to help achieve better bond strength after try-in*

Clinical Significance

The use of ZirClean will help to achieve reliable cementation results with Zirconia restorations after intraoral try-in.*



Shear Bond Strength to Zirconia[†] (MPa) After Cleaning With...



[†]Bonding after saliva contamination using Ultrabond shear bond test using Z-Prime™ Plus and Duo-Link Universal™

RELATED MATERIALS



ALL-BOND UNIVERSAL®



Z-Prime™ Plus

ORDER INFORMATION

1 Syringe (5g), 2 Disposable Syringe Tip Pack (15pk),
Instructions **B-7351P**

*As compared to untreated samples. Data on file.
Rx Only

Exclusively distributed by Curion

curion.ca | EN: 1.800.667.8811 | FR: 1.800.211.1200

 **curion**
formerly Bisco Canada