

LM-ARTE™



LM-Arte[™] instruments for aesthetic restorations

 $\mathsf{LM}\text{-}\mathsf{Arte}^{^\mathsf{TM}}$ is a range of innovative instruments for aesthetic restorations that is designed especially for composite layering. Each instrument is colour-coded and named after its main function, enabling you to administer easy, efficient and quick restorative treatments. The lightweight and ergonomic LM™ handle design, together with high-quality tip materials, provides excellent tactile sensitivity and easy handling for the best clinical results.

Reg. \$80.00

Please visit curion.ca to view current promotions.











LM-ARTE™ Dark Diamond





New LM-Arte[™] with LM Dark Diamond[™] coating

LM Dark Diamond[™] exclusive non-stick coated instruments are ideal for placing and contouring of composite resins. The dark color of the coating offers minimal reflections, and maximum contrast against composite. Hard and smooth surface allows for pleasant and precise restorative instrumentation protecting the tip from scratches and being easy to clean.

Reg. \$99.00


Please visit curion.ca to view current promotions.



LM Dark Diamond™ offers unique benefits!

- Contemporary diamond-like carbon (DLC) coating reduces stickiness for smooth composite placement
- Hard surface protects the tip from scratches
- Smooth surface for easy cleaning
- Dark coating color for minimal reflections, and maximal contrast against fillings
- Ergonomic LM-ErgoSense









LM Dark Diamond™



Non-Stick Coated Instruments

LM Dark Diamond instruments with non-stick coating for easy placing and contouring of composite resins. The diamond-like coating reduces the stickiness and allows pleasant and precise restorative instrumentation.

Reg. \$88.00

Please visit curion.ca to view current promotions.





Plastic Instrument - Mini Universal

LM444-445DDES

• For placing and contouring composite



Plastic Instrument - Anterior

LM43-47DDES

• For placing, packing and contouring composite



Plastic Instrument - Posterior

LM48-702DDES

• For placing, packing and contouring composite



Burnisher

LM31-32DDES

• For packing and modeling composite



Burnisher-Plugger

LM471-473DDES

• For modeling composite, Ball-ended, 1.5-2.5 mm



Packer-Modeller IO

LM482-483DDES

• For packing and modelling composite



Modeller MH

LM708-709DDES

• For packing large approximal fillings



Modeller MH Mini

LM707-708DDES

• For packing small and medium approximal fillings





LM Supportive Instruments





LM Arte Gengiva Kit

Obtain the best clinical results while avoiding unnecessary injury to the gingiva.

- Strong, yet flexible plastic tips adapt tightly to the gingival margin with minimal trauma to soft tissue
- Protects soft tissues from burs and abrasive rubbers
- Allows finishing instruments to reach the margins of the restorations
- Minimizes potential for debonding
- Easy to modify/mold edges to better fit the tooth to be restored
- Disposable

Kit includes:

- · LM-Gengiva Hand
- · LM-Gingival Retractors (blue)
- · LM-Multilever

Reg. \$138.00 Please visit curion.ca to view





GENGIVA

The gingival retracting kit





Retraction can be traumatic for the gingiva when not done properly. The use of LM-Gingival Retractors with LM-Gengiva instrument has many advantages over the traditional rubber dam and retraction cords.

LM-Gingival Retractors - Assortment of three sizes (90 pcs, 30 pcs/size): S for incisors, M for premolars and XL for molars.



Easy to customize

Edges can easily be molded in order to better fit to the tooth to be restored.



Better adaptation

The plastic tips allow enough flexibility for more adaptation while pressured.



Superior retraction and protection

The retractor is optimally rigid to maximize retraction while protecting the soft tissue from burs and abrasive rubbers



Improved reach to margins

Allows finishing instruments to reach the margins of the restoration while minimizing soft tissue trauma



Exclusively distributed by Curion