

Profiller Capsule Machine Comparison



Model	Profiller 1100	Profiller 3600	Profiller 3700	Profiller 3800
Capsules Filled	100	300	300	300
Handheld Base Model	Yes	Yes	No, Handheld Orienter	No, European Orienter
Capsule Sizes Available	00, 0, 1, 2, 3, 4, 5 (100 holes) 000 (80 holes)	00, 0, 1, 2, 3, 4, 5 (300 holes) 000 (180 holes)	00, 0, 1, 2, 3, 4, 5 (300 holes) 000 (180 holes)	00, 0, 1, 2, 3, 4 (300 holes)
Output	2000-3000 per hour	4600-6900 per hour	7200-9000 per hour	7200-9000 per hour
Types of capsules	Gelatin, HPMC (Veggie)	Gelatin, HPMC (Veggie)	Gelatin, HPMC (Veggie)	Gelatin, HPMC (Veggie)
Dimensions (W x D x H)	8 x 7 x 4 in	14.2 x 10.5 x 5 in	21 x 17.7 x 24 (filler) 11.5 x 7.8 x 2.4 (orienter)	21 x 17.7 x 24 (filler) 16.3 x 16.5 x 18.5 (orienter)
Cleaning	Dishwasher and automatic dishwasher detergent safe			
Chemicals Known to Discolour Parts	Nylon parts are known to discolour when used with Erythrosine. The powder spreader is known to discolour when used with Beta Carotene or Sunset Yellow	Nylon parts are known to discolour when used with Erythrosine. The powder spreader is known to discolour when used with Beta Carotene or Sunset Yellow	Nylon parts are known to discolour when used with Erythrosine. The powder spreader is known to discolour when used with Beta Carotene or Sunset Yellow	Nylon parts are known to discolour when used with Erythrosine. The powder spreader is known to discolour when used with Beta Carotene or Sunset Yellow
CE Certification	A declaration of conformity is included in each system	A declaration of conformity is included in each system	A declaration of conformity is included in each system	A declaration of conformity is included in each system
Shipping Weight	14 lbs	43 lbs	28 lbs and 103 lbs	29 lbs, 103 lbs and 116 lbs

Importers and Distributors of Compounding Products

PO Box 8009, Gold Coast Mail Centre
 QLD 9726, Australia

T: +61 (0)7 5597 4169
 F: +61 (0)7 5597 4264

E: bellacorp@bellacorp.com.au
 W: www.bellacorp.com.au