



Innovation and more!

An overview of the new developments of the gako unguator




 quick-lock coupling:
The mixing shaft will be loaded directly from the front into the lifting arm.




 Slide the mixing unit upwards...




 ...and with a slight turn the unit will be fixed in place - ready!




 The lifting arm can be easily and quickly replaced - e.g. to stir 2000 ml preparations.



 Connectivity to accommodate, electronic scales, printers, network, keyboard and more



 State of the art touch technology even for the BASIC and EMP for easy handling and cleaning.



Discover the quick-lock coupling following the link: goo.gl/gtDmLt