

5 REASONS TO TRY CORE-FLO™ DC CORE-FLO™ DC

Dual-Cure Core Build-Up, Endodontic Post Cement & Dentin Replacement Material

Rx Only

STACKABLE OR SELF-LEVELING



And low water sorption assures a strong foundation for crowns and direct restoratives.

2 Excellent Flowability

Void-free flow makes it ideal for post cementation.

3 Low Volumetric Shrinkage

Properties allow it to be used effectively as a dentin replacement material - cuts just like dentin.

4 Radiopaque

Formula makes it easy to identify on radiographs.

5 Versatility

With dual-cure chemistry that performs equally well in light-cure and self-cure modes.



MC-51210CF

WORKS WITH OR WITHOUT POSTS!



Exclusively distributed by Curion

curion.ca | EN: 1.800.667.8811 | FR: 1.800.211.1200