

## Choose sharp.

### Shape and Safety:

Sharpen-free instruments consistently maintain their original design, shape, and angulation throughout their entire lifespan. Excessive sharpening of instruments can alter the original cutting shape and angularity, potentially compromising instrument durability and increasing the risk of breakage.

### Savings Over Time:

Transitioning from conventional to sharpen-free instruments in a small clinic can lead to significant cost reductions, potentially saving up to \$10,000 per year. By opting for sharpen-free hand instruments, the need for sharpening is eliminated, saving both time and expenses, which can lead to an increase in productivity.








Opening a new office or upgrading multiple kits?  
Contact us for bulk pricing options.

Reg. \$100 - \$116


Please visit [curion.ca](http://curion.ca) to view current promotions.








### Sickle Scalers



### Gracey Curettes



### Universal Curettes





## Have confidence in your edge.

Experience the exceptionally ergonomic benefits of LM-DuraGradeMax's lasting sharpness, and durability when removing dental deposits. LM's combination of advanced edge retention technology and the highest-rated ergonomic handles on the market means less frequent sharpening and easier, more comfortable scaling for both the clinician and patient.

Opening a new office or upgrading multiple kits?  
Contact us for bulk pricing options.

Reg. \$76 - \$85

Please visit [curion.ca](http://curion.ca) to view current promotions.



Sickle Scalers		Dual Graces		Curettes		Gracey Curettes		Specialty	
	Nebraska 128 Sickle Scaler LM128-128ES		M23 Sickle Scaler LM313-314ES		Standard Dual-Gracey (Syntette) LM215-216ES		Mini Gracey 7/8 Curette LM207-208MES		Standard Gracey 15/16 Curette LM245-246ES
	H6/H7 Sickle Scaler LM146-147ES		204S Sickle Scaler LM303-304ES		Mini Dual-Gracey (Syntette) LM215-216MES		Standard Gracey 11/12 Curette LM211-212ES		Mini Gracey 15/16 Curette LM245-246MES
	Laplander Sickle Scaler LM148-149ES		Barnhart 5/6 Curette LM227-228ES		Standard Gracey 1/2 Curette LM201-202ES		Mini Gracey 11/12 Curette LM211-212MES		Mini Gracey 17/18 Curette LM247-248MES
	Micro Sickle Scaler LM301-302ES		Columbia 4L-4R Curette LM219-220ES		Mini Gracey 1/2 Curette LM201-202MES		Standard Gracey 13/14 Curette LM213-214ES		Standard Gracey 17/18 Curette LM247-248ES
	Mini Sickle Scaler LM311-312ES		McCall 13S-14S Curette LM222-223ES		Standard Gracey 7/8 Curette LM207-208ES		Mini Gracey 13/14 Curette LM213-214MES		Furcator KS LM259-260ES

## Comfortable and lightweight.

Each LM diagnostic instrument features the highest-rated ergonomic handle in the industry. Made of durable medical-grade silicone, LM-Instruments are comfortable, lightweight, and won't fade or change with repeated cleaning cycles. Every tip is made of LM's proprietary wear-resistant DuraGradeMax supersteel.

LM ErgoSense

Opening a new office or upgrading multiple kits?  
Contact us for bulk pricing options.

Reg. \$18 - \$65

Please visit [curion.ca](http://curion.ca) to view current promotions.



### HYGIENISTS' TOP PICK



#### Flexplorer 11F-12F

LM11F-12FES

- Flexible calculus detection explorer angulated similarly to Gracey 11/12
- Well-suited for periodontal diagnosing of subgingival calculus and probing deep pockets and furcations



#### Explorer - Periodontal Probe 23-520B

LM23-520BES

- Combination instruments for probing caries and parodontium.
- 2 mm scale



#### Marquis 12 Probe 53B

LM53BES

- Instrument for probing caries and parodontium.
- 3 mm scale



#### Nabers Probe

LM20-2IES

- For probing the extent and depth of furcation lesions.



#### Mirror Handle 28

LM28ESGREY

- Cone socket LM 28 mirror handle threads are in accordance with the international standard UNF4-48.

### DOCTORS' TOP PICK



#### Explorer Centered 17C-23

LM17C-23ES

- For diagnosing calculus and caries and exploring of pockets, restorations and furcations.
- Especially for probing distal surfaces



#### Explorer - Periodontal Probe 23-530B

LM54BES

- Combination instruments for probing caries and parodontium.
- 3 mm scale



#### Furcation Probe

LM20B-21BES

- For probing the extent and depth of furcation lesions.
- 2 mm scale



#### North Carolina Probe 54B

LM54BES

- Combination instruments for probing caries and parodontium.
- Scale 5-10-15



#### Williams Probe

LM5IES

- Instrument for probing caries and parodontium.
- Scale 1-2-3-5-7-8-9-10



## Safety and functionality.

State of the art series of titanium implant maintenance instruments expertly designed to ensure excellent tactile sensitivity while gently, but effectively removing bacterial biofilm and calculus deposits on implant surfaces without harming abutments.

Opening a new office or upgrading multiple kits? Contact us for bulk pricing options.

**\$572**  
Kit

**\$135**  
Instrument

**\$34**  
Refill Tip

Please visit [curion.ca](http://curion.ca) to view current promotions.



### A great instrument series for implant maintenance.

The LM-ErgoMix™ structure ensures excellent tactile sensitivity and optimal rigidity. Titanium tips for sensitive but still effective implant maintenance.

### Benefits of retipping:

- Always sharp instruments
- Minimizes waste
- Enables personalized tip combinations



#### Implant Mini Universal Curette

LM283-284MTIEM

- Delicate working end
- For premolars and molars in the upper jaw
- Especially for tight interdental spaces.
- Better access into tight deep pockets and furcations

**Refill Tip**  
Right Tip: LM284MTIEM  
Left Tip: LM283MTIEM



#### Implant Mini Gracey 1/2

LM201-202MTIEM

- For incisors and premolars

**Refill Tip**  
Mini Gracey 1 Tip: LM201MTIEM  
Mini Gracey 2 Tip: LM202MTIEM



#### Implant Mini Gracey 11/12

LM211-212MTIEM

- For premolars' and molars' mesial surfaces

**Refill Tip**  
Mini Gracey 11 Tip: LM211MTIEM  
Mini Gracey 12 Tip: LM212MTIEM



#### Implant Mini Gracey 13/14

LM213-214MTIEM

- For premolars' and molars' distal surfaces

**Refill Tip**  
Mini Gracey 13 Tip: LM213MTIEM  
Mini Gracey 14 Tip: LM214MTIEM