



# GLUE RENEW

AN EASY WAY TO RECONDITION CLIMBING SKINS

- Reconditions old adhesives
- Suitable for all G3 skins
- G3 factory formulation
- Ultra-low temperature tolerant adhesive
- Manageable at very low temperatures

## CONTENTS OF GLUE RENEW PACKAGE:

2 lengths of glue renew product  
Re-gluing instructions

G3 Genuine Guide Gear Inc.  
[www.genuineguidegear.com](http://www.genuineguidegear.com)





Scan the QR code for video instructions, product information and more.

# GLUE RENEW

AN EASY WAY TO RECONDITION CLIMBING SKINS

## EQUIPMENT YOU WILL NEED:

- Working bench
- Waxing iron or household iron (no steam) - adjustable to 180°C
- Metal Paint Scraper
- Newspaper or cardboard
- Utility knife
- Clamps (2) (to hold skin to work bench)

## DIRECTIONS FOR REMOVING OLD ADHESIVE

Thoroughly remove old or contaminated climbing skin adhesive before applying G3 Glue Renew. G3 recommends the following method:

1. Lay skin adhesive side up on the workbench and clamp skin taut at both ends
2. If you have a heat gun, remove the old adhesive by moderately heating and scraping the adhesive onto a piece of cardboard. Skip to step 5. If you don't have a heat gun, cover the skin adhesive with 30cm long sections of old newspaper
3. Heat newspaper sections of the adhesive (approx. 30 cm at a time) at ~160°C by slowly passing the waxing iron over the newspaper
4. Immediately after heating, scrape the adhesive and newspaper off the skin with your metal paint scraper
5. Repeat process until no adhesive remains on the surface of the skin backing, discard old adhesive and used newspaper

## DIRECTIONS FOR APPLYING THE NEW ADHESIVE:


1. Peel back the Glue Renew plastic liner (transparent side) to expose the adhesive. Place the Glue Renew adhesive side down on the freshly scraped skin. Use a utility knife to cut to the desired length and width, take care not to cut the skin underneath
2. Iron Glue Renew onto the skin fabric at medium temp (160°C). Work from tip to tail applying medium pressure and moving slowly down the skin to evenly apply the adhesive. The pressure should be enough to even out the surface without the adhesive seeping out the sides. Repeat motion until paper backing appears smoother and more saturated. The adhesive may have small wrinkles and air bubbles, work in sections applying higher pressure for 2-3 seconds to even out these areas
3. At high temp (180°C), quickly pass over the adhesive surface and perimeter with light-medium pressure
4. Trim any excess adhesive using a knife and allow the adhesive to settle (G3 recommends 16-24hrs) before removing the release paper
5. After the adhesive has settled, peel up about an inch of paper at the skin tip in one continuous piece
6. Flip the skin so the plush is facing up and hold edge of the paper liner down with your palm while slowly pulling the skin away from the paper. If the paper rips, slowly peel at a 45° angle until it is a continuous piece again

### G3 PRO TIPS:

If the Glue Renew isn't fully embedded into the backing, cover the area with a piece of parchment paper and re-apply the Glue Renew at 180°C with high pressure. Allow the area to cool to room temperature and remove the parchment paper.

Excess glue can be cleaned off tools and surfaces by heating and scraping, using citrus base cleaner in a well-ventilated area, or by using an extra piece of Glue Renew to "pick up" glue.

Remaining G3 Glue Renew can be used as touch up adhesive, if required.



DESIGNED BY G3 IN CANADA  
CONÇU PAR G3 AU CANADA

MADE IN USA  
FABRIQUÉ AU USA

(418-620)  
G3-SKN-G005-REV-D