




ThermoWorks

# WAND™

Non-Contact, Digital Forehead Thermometer



**Non-contact, 1-second readings for quick, safe screening.**

*\*FDA-cleared for safety & accuracy for kids and adults. ThermoWorks Wand Thermometer measures infrared energy radiated from the skin at the center of the forehead area. This captured energy is collected through the lens and converted to a body temperature value (displays oral equivalent temperature).*

*It's helpful to know each individual's normal temperature when they are well. This is the best way to accurately diagnose a fever. To get a reference value, we suggest taking multiple readings (5 to 10) over one minute, taking the highest number to determine "normal temperature".*