



2000 ml / 2600 ml = Nominal-Volume / Filling-Volume

measurements in mm

Color Versions:

gako unguator jar w / white lid 2000ml AD

PCN:

AT1000000275

PPN:

111391832018

GTIN:

4251208804990

Material:

Product	Material	Colour concentration	
		standard	cosmetic
gako unguator Jar Cap	Polypropylene natural	white	white
gako unguator Jar Lid	Polypropylene natural	red, white, green, blue	white
gako unguator Jar Body	Polypropylene natural	white	white
gako unguator Jar Bottom	Polyethylene	-	-

The following five tests were carried out according to the guideline of the Federal Chamber of Pharmacies „testing and storage of primary packaging material“:

- Microbiological contamination (carried out according to Ph. Eur 2.6.12 and 2.6.13.)
 - Total number of colonies CFU \leq 1 per ml nominal volume
 - Exclusion of pathogenic indicator microorganisms
- Particle contamination (carried out according to DAC test 17)
 - All size classes: \leq 1000
 - 150 to 400 μm : \leq 100
 - $>$ 400 μm : \leq 100
- Extractable substances (color fastness) (carried out according to JP 7.02)¹
 - No change in absorbance between 220 and 350 nm
- Opacity (carried out according to Ph. Eur 3.2.1.)¹
 - Transmittance: $<$ 1% (in the spectral region from 290 to 450nm)
- Water vapor permeability (seal integrity) (carried out according to USP 671)
 - Water vapor permeability: \leq 100 mg per day and 1l nominal volume

¹ Not for transparent jars!

gako international

80687 Munich / Germany

Matthias Konietzko

Dipl.-Wirtsch.-Ing. Matthias Konietzko

2018-03-01