

Ordering

500021 – 60 mL

500886 – 500g

Applications

Suspensions

Product Highlights

SyrSpend® SF Alka dry is an alkaline suspending base for reconstitution to compound elegant oral liquid dosage forms with a pleasant taste. Available in validated, graduated dispensing containers allowing for an easy three-step compounding process directly into the dispensing container. This unique vehicle uses an active suspending technology to offer more accurate dosing by stabilising the suspension at rest and actively preventing 'caking' of insoluble ingredients at the bottom of the bottle.

- The end product is a ready-to-use, all-in-one suspending and sweetening liquid vehicle
- Compatible with an extensive range of APIs
- Buffered to pH >7.0 for maximum stability of acid-labile APIs
- Low osmolality (<50 mOsmol/kg)
- Preservative-free

Ingredients

Calcium Carbonate, Corn Starch, Sucralose, Cherry Flavor

Specifications

- pH: >7.0 (when reconstituted)

Formulas

#FA01122 Esomeprazole 3 mg/mL Oral Suspension in SyrSpend SF ALKA Dry

#FA00937 Furosemide 10 mg/mL Oral Suspension in SyrSpend SF ALKA Dry

#FA01017 Pantoprazole 2 mg/mL Oral Suspension in SyrSpend SF ALKA Dry

Formulas are available for download for FACTS members: fagronacademy.us